



**QUEEN'S
UNIVERSITY
BELFAST**

Loop Mediated Isothermal Amplification (LAMP) assay for Rapid detection of *Streptococcus agalactiae* (Group B *Streptococcus* - GBS) in vaginal swabs - A Proof of Concept Study

McKenna, J., Cox, C., Fairley, D. J., Burke, R., Shields, M. D., Watt, A., & Coyle, P. V. (2017). Loop Mediated Isothermal Amplification (LAMP) assay for Rapid detection of *Streptococcus agalactiae* (Group B *Streptococcus* - GBS) in vaginal swabs - A Proof of Concept Study. *Journal of Medical Microbiology*, 66(3), 294-300.
<https://doi.org/10.1099/jmm.0.000437>

Published in:
Journal of Medical Microbiology

Document Version:
Peer reviewed version

Queen's University Belfast - Research Portal:
[Link to publication record in Queen's University Belfast Research Portal](#)

Publisher rights
Copyright 2017 the authors.
This work is made available online in accordance with the publisher's policies. Please refer to any applicable terms of use of the publisher.

General rights
Copyright for the publications made accessible via the Queen's University Belfast Research Portal is retained by the author(s) and / or other copyright owners and it is a condition of accessing these publications that users recognise and abide by the legal requirements associated with these rights.

Take down policy
The Research Portal is Queen's institutional repository that provides access to Queen's research output. Every effort has been made to ensure that content in the Research Portal does not infringe any person's rights, or applicable UK laws. If you discover content in the Research Portal that you believe breaches copyright or violates any law, please contact openaccess@qub.ac.uk.

Loop-Mediated Isothermal Amplification (LAMP) assay for Rapid detection of *Streptococcus agalactiae* (Group B Streptococcus - GBS) in vaginal swabs – A Proof of Concept Study

J.P McKenna^{1*}, C. Cox¹, D.J. Fairley¹, R. Burke^{1,2}, M.D. Shields^{2,3}, A. Watt¹, P.V. Coyle¹

¹ Regional Virus Laboratory, Department of Microbiology, Belfast Health and Social Care Trust, Belfast, UK.

² Centre for Experimental Medicine, Queens University Belfast, UK.

³ Royal Belfast Hospital for Sick Children, Belfast Health and Social Care Trust, UK.

* Corresponding author. Email: james.mckenna@belfasttrust.hscni.net. Tel. +44(0)2890 635303

Abstract

Neonatal infection with *Streptococcus agalactiae* (Group B streptococcus – GBS) is a life threatening condition which is preventable if colonised mothers are identified and given antibiotic prophylaxis during labour. Conventional culture is time consuming and unreliable, and many available non-culture diagnostics are too complex to implement routinely at point of care. Loop-mediated isothermal amplification (LAMP) is a method which enables the rapid and specific detection of target nucleic acid sequences in clinical material without any requirements for extensive sample preparation. A prototype LAMP assay targeting the GBS *sip* gene is described. The assay was 100% specific for GBS with a Limit of Detection of 14 genome copies per reaction. The clinical utility of the prototype LAMP assay for rapid direct molecular detection of GBS was determined by testing a total of 157 vaginal swabs with

minimal sample processing using a rapid lysis solution. Compared to a reference qPCR assay the direct LAMP protocol had a sensitivity and specificity of 95.4% and 100% with a Positive Predictive Value (PPV) and Negative Predictive Value (NPV) of 100% and 98.3%. Positive (LR+) and negative (LR-) likelihood ratios were infinity and 0.05 respectively. The direct LAMP method took a mean of 45 minutes from receipt of a swab to generation of confirmed result compared to 2 hours 30 minutes for the reference qPCR test. The direct LAMP protocol described is easy to perform facilitating rapid and accurate detection of GBS in vaginal swabs. This test has potential for use at point of care.

Introduction

Streptococcus agalactiae (Group B streptococcus - GBS) is the leading cause of early onset neonatal sepsis [1-5]. Asymptomatic GBS colonisation of the vagina and lower gastrointestinal tract occurs in 10-30% of pregnant women and is the most significant risk factor for neonatal GBS disease as transmission can occur from the recto/vaginal area during childbirth [4,6-9]. The implementation of culture based antepartum screening programmes for all pregnant women at 34-38 weeks and initiating intrapartum antibiotic prophylaxis has helped to reduce but not eliminate neonatal early onset disease (EOD) in many industrialised countries [4,10,11]. Despite technological advances the primary method for detection of GBS colonisation remains selective laboratory based culture of vaginal and/or rectal swab specimens followed by ancillary confirmation tests. This approach demonstrates variable sensitivity [12-14] and is associated with long turnaround times of 18-72 hours rendering culture unfeasible as an intrapartum test. Thus, women who give birth prematurely are not screened, and women who become colonised after screening are not identified [7,8]. A sensitive screening assay which could rapidly determine GBS status during labour would be

beneficial as it would facilitate targeted therapy whilst reducing the risks associated with antibiotic use in non-colonised women [15].

Quantitative Polymerase Chain Reaction (qPCR) assays can provide test results in a few hours [15-20]. However, requirements for expensive equipment, intricate sample preparation, skilled staff and/or support precludes the use of PCR at point of care (POC) [21-23].

Loop-mediated isothermal amplification (LAMP) is an isothermal molecular amplification technique which offers simple qualitative detection combined with exceptional sensitivity and specificity [24-26]. The single reaction temperature (60 – 65°C) employed by LAMP means there are no requirements for thermocycling resulting in a need for less complex equipment making it amenable to use at POC [22,27]. Furthermore compared to PCR, LAMP is also demonstrably more resistant to inhibitors found in clinical material indicating potential for direct application to clinical specimens and/or crude sample preparations [21,28,29]. This is advantageous as one of the issues frequently hampering application of Nucleic acid amplification tests (NAAT) at POC are requirements to perform dedicated multistep nucleic acid extraction steps which add to test complexity, cost and turnaround time [21-23]. Whilst a number of LAMP assays for GBS detection have been developed previously, they are not amenable for POC use having requirements for an 18-24 hr pre-enrichment culture step [20], require multiple time consuming manipulations prone to contamination [30] or have unproven clinical efficacy [31,32].

Here we describe a newly designed GBS LAMP assay targeting *S. agalactiae sip* gene. As a proof of concept the LAMP assay was applied to residual vaginal swabs which had undergone minimal sample processing (Direct detection – no extraction protocol) with a rapid lysis solution containing the bacteriolytic enzyme achromopeptidase (ACH). ACH demonstrates high lytic activity against a range of gram positive cocci including GBS [23,33-35].

Materials and Methods

Clinical Specimens: The study was approved by the Office for Research Ethics committee Northern Ireland (reference 09/NIR02/43). A total of 157 residual vaginal swabs from pregnant women submitted to Department of Microbiology, Belfast Health and Social Care Trust (BHSCT) for bacterial investigations between March – July 2013 were opportunistically collected and utilised for the clinical validation study. The sample set consisted of 153 High Vaginal swabs (HVS) and 4 Low Vaginal Swabs (LVS). All vaginal swabs were collected on plastic rayon tipped swabs (Sterilin Ltd., Newport, UK) and transported to laboratory in AMIES transport medium. Culture of all swab specimens was carried out according to BHSCT Standard Operating Procedure by laboratory staff unaware that the samples were to be subsequently included in a study. Briefly, each swab was plated onto Columbia Blood agar (CBA) and differential Granada™ agar (bioMérieux, Basingstoke, UK). CBA plates were incubated in 10% CO₂ atmosphere with Granada™ plates incubated anaerobically at 37°C. Plates were examined at 18-24 hours and again at 48 hours with presumptive positive GBS colonies on Granada™ agar identified by production of orange pigment. Presumptive GBS colonies on CBA were identified on the basis of haemolysis. All suspect colonies were appropriately identified by a combination of standard tests including Gram stain, catalase activity and latex agglutination using SLIDEX® Strepto Plus kit (bioMérieux, Basingstoke, UK) to identify the presence of Group B antigen. Following routine bacteriological culture all residual swabs were stored at +4°C for a maximum of one week until required for processing and molecular analysis.

Rapid Sample preparation method for direct LAMP: Lysis solution consisted of 10 U/L of achromopeptidase (ACH) (Sigma-Aldrich, Dorset, UK) in 10 mmol/L Tris-HCL pH8.0, 1

mmol/L EDTA (Sigma-Aldrich, Dorset, UK). Residual swab specimens were removed from AMIES transport medium tubes with ~2 cm end physically broken off into 400 µl pre aliquoted lysis solution [ACH] by firmly pressing the swab against the side of a microfuge tube (Sarstedt, Leicester, UK). Tubes were capped and swabs incubated at ambient 'room temperature' (22 - 24°C) for 5 minutes before boiling at 95°C for 5 minutes in a dedicated heat block to simultaneously deactivate ACH and denature any target DNA for LAMP analysis. Samples were immediately cooled by placing directly on ice. A 5µl aliquot of lysis supernatant was used for GBS direct LAMP testing as outlined below. Remaining lysis supernatant was stored at 4°C.

Nucleic acid (NA) extraction for qPCR and indirect LAMP: GBS reference and clinical isolates were pre-treated with ACH before undergoing total nucleic acid extraction using QIAamp blood DNA mini extraction kit (Qiagen, Crawley, UK). A 200 µl saline suspension of each GBS isolate was prepared and added to an equal volume of ACH 20 U/L for a final working concentration of ACH of 10 U/L. The suspension was incubated at room temperature for 5 minutes before being boiled at 95°C/5 minutes. 200 µl aliquot of each lysate was then added to AL lysis buffer (Qiagen, Crawley, UK) and processed as per manufacturer's instructions to produce 100 µl of genomic extract.

Bacterial and fungal genomic DNA for specificity evaluation was obtained using QIAamp blood DNA mini extraction kit (Qiagen, Crawley, UK). 200 µl saline suspensions of each target organism were added to AL lysis buffer and processed as per manufacturer's instructions to produce 100 µl of genomic extract.

Vaginal swab lysates were extracted as follows; A 200 µl aliquot of ACH lysis supernatant (produced as indicated in previous section) was removed from each clinical swab specimen and extracted using MagNA Pure 96 viral NA small volume kit and MagNA Pure instrument

(Roche Diagnostics, Mannheim, Germany) as per manufacturer's instructions. Resultant 100 µl genomic extract aliquots were stored at -80°C until required for qPCR and LAMP (indirect) analysis.

LAMP primers and reaction: LAMP primers targeting a highly conserved region of *S. agalactiae sip* gene (HiberGene Diagnostics, Ireland) were employed throughout. The LAMP primer set consisted of forward (sip_FIP) and reverse (sip_BIP) 'inner' primers, forward (sip_F3) and reverse (sip_B3) 'outer' primers and forward (sip_FL) and reverse (sip_BL) 'loop' primers. Optimised GBS LAMP consisted of 5 µl of template DNA, 12.5 µl of LAMP isothermal Master Mix ISO 001 (OptiGene Ltd., Horsham, UK) containing novel strand displacing DNA polymerase (*GspSSD*), propriety thermostable inorganic pyrophosphatase, reaction buffer, MgCl, dNTP's, DNA intercalating dye, 2 µM each of FIP and BIP, 0.2 µM each of F3 and B3, 1 µM each of FL and BL and Nuclease Free Water (Promega, Southampton, UK) for a final reaction volume of 25 µl. The GENIE® II (OptiGene Ltd, Horsham, UK) real time fluorometer was utilised for concurrent target isothermal amplification and detection. LAMP reactions were performed at a previously identified optimal reaction temperature of 62°C for 30 minutes followed by a heating and cooling step to 98–80 °C (0.05 °C/s) which allows re-annealing of amplified DNA and display of a specific annealing curve. Positive and negative controls were included for each run and stringent precautions undertaken to prevent cross contamination. Positive LAMP reactions were identified by GENIE® II (OptiGene Ltd, Horsham, UK) instrument software. In the results section each positive reaction had a Time to positivity (T_T) value in minutes and seconds and a specific melting point temperature (T_m) in degrees centigrade (°C).

Real Time PCR (qPCR): qPCR was employed as the reference molecular detection method. Primer and hybridisation probe sequences targeting *S. agalactiae sip* gene were as previously described [17]. Primer sequences as follows: Forward primer 5'–ATCCTGAGACAACACTGACA–3', and reverse primer 5'–TTGCTGGTGTCTTCTATTTTCA–3'. The hybridisation probe 5'–ATCAGAAGAGTCATACTGCCACTTC3' contained a FAM label on 5' end with a Black Hole Quencher 1 (BHQ1) on 3' end. Primers and hybridisation probe were synthesised by Eurogentec (Eurogentec, Liège Belgium). The final qPCR reaction mix contained 1X Platinum® UDG Mastermix (Life Sciences, Paisley, UK), 0.2 µM Bovine serum albumin (Sigma, Dorset, UK), total of 4 mmol/L MgCl₂, 0.4 µM forward and reverse primers, 0.2 µM hybridisation probe, Nuclease Free Water (Promega, Southampton, UK) and 3 µl of target template for a final reaction volume of 12 µl. qPCR was performed using a Light Cycler 480 (LC480) instrument (Roche Diagnostics, Mannheim, Germany). A No Template Control (NTC) consisting of Nuclease Free Water (Promega, Southampton, UK) and a plasmid reference standard containing 3.50E06 *sip* genome copies per ml were included for all qPCR runs. The qPCR run parameters consisted of an initial step of 95°C for 10 minutes followed by an amplification program for 45 cycles of 10 seconds at 95°C, 30 seconds at 60°C, 1 second at 72°C with fluorescence acquisition at the end of each extension. Data was analysed using LC480 software and GBS genome copy number for positive specimens determined from crossing point threshold (Cp) relative to an external positive control curve.

Specificity & sensitivity: Specificity of GBS LAMP was evaluated using total genomic extracts (containing approximately 10E07 genome copies) derived from a range of bacteria (n=55) and fungi (n=1). These included a number of GBS reference, clinical and veterinary isolates (n=17). Reaction mix without DNA template was included as a negative control.

174 The following organisms were tested; *Streptococcus agalactiae* (NCTC 11360, NCTC 8187,
 175 NCTC 11078, NCTC 11079, NCTC 11080, NCTC 11930, NCTC 9829, ATCC 12386, ID:
 176 603532, ID: 603535, ID: 605872, RVH2419/12, RVH2543/12, RVH60/13, RVH524/13,
 177 VSD04988, VSD04946), *Streptococcus pneumoniae* (ATCC 49619, RVH59/07),
 178 *Streptococcus parasanquis* (RVH74/07), *Streptococcus intermedius* (RVH305/07),
 179 *Streptococcus pyogenes* (ATCC 19615), *Streptococcus salivarius* (RVH557/13),
 180 *Streptococcus vestibularis* (RVH527/13), *Streptococcus mitis/oralis* (RVH508/13),
 181 *Streptococcus dysgalactiae* (RVH497/13), *Streptococcus anginosus* (RVH452/13),
 182 *Staphylococcus aureus* (ATCC 25923), MRSA (NCTC 10442), *Staphylococcus epidermidis*
 183 (ATCC 12228), *Staphylococcus lugdunensis* (RVH isolate), *Staphylococcus capitis* (RVH
 184 isolate), *Enterococcus faecalis* (ATCC 29212, RVH009/07), *Micrococcus luteus*
 185 (RVH290/07), *Brevibacillus brevis* (RVH381/13), *Bacillus cereus* (RVH299/13), *Neisseria*
 186 *meningitidis* (NCTC 10791 serogroup Y, RVH531/07 serogroup B), *Neisseria sicca*
 187 (RVH920/06), *Klebsiella pneumoniae* (RVH45/07), *Klebsiella oxytoca* (RVH135/07),
 188 *Enterobacter aerogenes* (NCTC 10006, RVH001/07), *Enterobacter cloacae* (NCTC 9394),
 189 *Serratia marcescens* (RVH287/07), *Acinetobacter baumannii* (RVH36/07), *Acinetobacter*
 190 *lwoffii* (RVH98/07), *Escherichia coli* (NCTC 9001, ATCC 25922), *Citrobacter freundii*
 191 (NCTC 9750), *Haemophilus influenzae* (NCTC 4560), *Sphingobacterium spiritivorum*
 192 (ATCC 33861), *Stenotrophomonas maltophilia* (ATCC 17666), *Pseudomonas aeruginosa*
 193 (ATCC 27853), *Gardnerella vaginalis* (ATCC 49145), *Candida albicans* (ATCC 10028).
 194 LAMP analytical sensitivity (limit of detection - LOD) was determined using tenfold serial
 195 dilutions of *S. agalactiae* ATCC 12386 nucleic acid extracts analysed in duplicate containing
 196 2.8E07 to 2.8E02 genome copies per ml respectively. A No Template Control (NTC)
 197 consisting of Nuclease Free Water (Promega, Southampton, UK) was included with all
 198 LAMP runs.

Statistical analysis: For purposes of validation the reference molecular method was defined as achromopeptidase (ACH) pre-treatment of a clinical swab specimen followed by Nucleic Acid (NA) extraction and qPCR testing. A specimen was considered a true positive if: (i) GBS was obtained following culture; or (ii) Specimen was identified as GBS qPCR positive. All qPCR and LAMP results were compared with routine culture results. Sensitivity, specificity, positive predictive value (PPV), negative predictive values (NPV), positive (LR+) and negative (LR-) likelihood ratios for each test method were calculated. The level of agreement between LAMP and qPCR assays was assessed using Kappa (κ) coefficient with 95% confidence levels. κ was calculated using the formula $\frac{\text{Pr}(o) - \text{Pr}(e)}{1 - \text{Pr}(e)}$, where $\text{Pr}(o)$ is the observed agreement and $\text{Pr}(e)$ is the hypothetical probability of chance agreement.

Results

Analytical reactivity, specificity and sensitivity: The LAMP assay successfully amplified DNA from all GBS reference, clinical and veterinary strains analysed with Time to positivity (T_T) values of less than 20 minutes and specific melt point temperature values (T_m) ranging from 85.5 – 86.21⁰C obtained. No amplification was recorded for any other bacterial or fungal species tested. This confirmed 100% specificity of the GBS LAMP assay.

Employing tenfold serial dilutions of a known concentration (2.8E07 to 2.8E02 copies per ml) of GBS ATCC 12386 nucleic acid extracts in duplicate the LAMP assay consistently detected 2.80E03 *sip* genome copies per ml corresponding to an analytical Limit of Detection (LOD) of 14 genome copies per reaction (for 5 μ l input volume and final reaction volume of 25 μ l). The mean GBS LAMP T_T and T_m for 2.80E03 *sip* genome copies per ml were determined as 11 minutes and 85.98⁰C respectively. The LAMP LOD was slightly higher than that achieved for qPCR which in our laboratory had a previously determined LOD of 1.00E03 genome copies per ml equating to 2 genome copies per reaction.

Clinical Validation: Routine bacteriological culture was compared to the reference qPCR molecular method, direct LAMP and indirect LAMP methods. Summary of findings are outlined in Tables 1 and 2. A total of 36/157 (22.9%) specimens analysed were confirmed as GBS culture positive. All 36 culture positive GBS specimens were correctly identified as positive by qPCR and by the direct and indirect LAMP methods. An additional 7 specimens were identified as GBS positive by qPCR giving a total of 43/157 (27.4%) ‘true’ positive specimens. Comparison of difference between qPCR and culture results was statistically significant ($p=0.0005$). The loads of GBS determined by qPCR (genome copies per ml) for positive clinical specimens were as follows: Mean = 1.23E+08; Max= 2.50E+09; Min=

249 1.90E+03; Median= 2.40E+06. Compared to culture qPCR had a sensitivity and specificity of
 250 100% and 94.2% with a PPV and NPV of 83.7% and 100%. Positive (LR+) and negative
 251 (LR-) likelihood ratios were infinity and 0.16 respectively. The mean genome copies per ml
 252 of the seven identified specimens missed by culture but positive by qPCR was 1.09E06
 253 (Range 1.90E03 – 7.60E06; Median 5.50E03) compared to 1.47E08 (Range 4.60E03 -
 254 2.50E09; Median 3.65E06) for specimens identified as both culture and qPCR positive. This
 255 difference in mean genome copies per ml was statistically significant ($p < 0.0001$).
 256 The direct LAMP method took a mean of 45 minutes from receipt of a swab to generation of
 257 confirmed result compared to 1 hour 45 minutes for indirect LAMP (extraction protocol) and
 258 2 hours 30 minutes for reference qPCR test.
 259 There was a very high level of agreement between LAMP (Indirect and direct) methods and
 260 reference qPCR assay. Compared to qPCR agreement ranged from 100% for indirect LAMP
 261 (Kappa = 1.0; Strength of agreement = perfect) to a level of agreement of 98.7% for direct
 262 LAMP (Kappa = 0.968; Strength of agreement = very good). Indirect LAMP method
 263 identified 43/157 (27.4%) specimen extracts as GBS positive with the direct LAMP method
 264 identifying a total of 41/157 (26.1%) of crude specimen ACH lysates as GBS positive.
 265 Compared to qPCR the indirect LAMP method had a sensitivity and specificity of 100% and
 266 a PPV and NPV of 100%. Positive (LR+) and negative (LR-) likelihood ratios were infinity
 267 and 0.00 respectively. Compared to qPCR the direct LAMP had a sensitivity and specificity
 268 of 95.4% and 100% with a PPV and NPV of 100% and 98.3%. Positive (LR+) and negative
 269 (LR-) likelihood ratios were infinity and 0.05 respectively. Two specimens (405 & 453 see
 270 Table 2) were identified as false negative after analysis by direct LAMP method. Following
 271 nucleic acid extraction and analysis by qPCR these specimens tested GBS positive with
 272 genome copy loads of 1.90E+03 per ml (9.5 copies per 5µl addition) and 1.91E03 per ml
 273 (9.6 copies per 5µl addition) respectively. Additionally following nucleic acid extraction both

of the specimens tested LAMP positive returning T_T values of 12.15 and 29.00 minutes respectively. Overall the mean time to achieve positivity (T_T) for direct LAMP method was higher compared to indirect LAMP method; 14.83 minutes (range 8.15 – 26.00 minutes) compared to 9.60 minutes (range 5.00 – 29.00 minutes). Testing performance comparison is displayed in Table 1 with all specimens identified as positive by either culture, qPCR and/or LAMP is presented in Table 2.

Discussion

Screening for maternal GBS colonisation at 35-37 weeks via selective culture is a time consuming and inadequate method which fails to identify all cases of GBS [12,13,20]. There is a recognised clinical need for rapid, highly sensitive and specific methods to detect GBS which can be used at POC before or during labour to identify women for prophylactic antibiotic treatment. Currently qPCR offers the required specification to achieve these goals however technical obstacles including requirements for stringent sample preparation, complex instrumentation and skilled personnel have restricted their introduction at POC [21-23]. The robust chemistry, ease of operation, rapidity, low cost, and modest equipment requirements associated with LAMP represent significant advantages over existing molecular methods for detection of pathogens [22,27,29,36].

The primary objective of this proof of concept study was to validate the effectiveness of a prototype LAMP assay (HiberGene diagnostics, Ireland) in combination with a simple rapid specimen lysis protocol (without extraction) for direct detection of GBS in residual vaginal swabs. Pre-lysis with ACH at a concentration of 10 U/L for 5 minutes at room temperature has been previously shown in our laboratory to be a highly effective method to significantly increase the available yield of GBS nucleic acid for molecular analysis [35]. Eliminating the

requirement for nucleic acid extraction not only simplified the test procedure but also resulted in significant time savings with direct LAMP taking 45 minutes from receipt of swab to definitive GBS result compared to 60 minutes and >2 hrs for indirect LAMP and qPCR respectively. The GBS LAMP assay had an analytical LOD of 2.80E03 genome copies per ml which compared favourably with the determined reference qPCR analytical LOD of 1.00E03 genome copies per ml. Crucially when applied to clinical specimens, qPCR and both direct and indirect LAMP methods resulted in identification of additional GBS positives (n=7) missed by culture. These findings were not unexpected and are consistent with previous studies demonstrating an increase in GBS positivity rate when Nucleic acid amplification tests (NAAT's) are directly compared to culture [13,16,18,19,20]. Higher detection rates of NAAT's have been claimed to be due to the inability of culture to detect low numbers of GBS bacteria but may also be attributed to the presence of antagonistic flora and/or the detection of non-culturable cells [19,37].

Compared to the reference qPCR method the direct LAMP had a sensitivity of 95.4% and specificity of 100%. High specificity is an intrinsic feature of LAMP and reflects the presence of six specific LAMP primers recognising eight target regions on *sip* gene. The direct LAMP assay returned higher mean T_T values for identified positive specimens compared to the indirect LAMP assay applied to the same samples which had undergone nucleic acid extraction (14.69 minutes versus 9.61 minutes respectively - difference of 5.08 minutes). Additionally the direct LAMP assay failed to detect GBS in two specimens subsequently shown to contain low genomic loads by qPCR. These findings indicate that the presence of component(s) in crude lysate matrices can impede LAMP and in conjunction with low initial GBS loads may have contributed to the two false negative results. Solution(s) to combat this could include prior centrifugation and/or a dilution step(s) although this would complicate the assay protocol and make application less attractive at POC. At present the

clinical significance of specimens containing a low GBS burden remains to be determined [16]. The excellent performance of qPCR and LAMP in this study must also be taken in context as ‘second hand’ residual vaginal swabs post bacteriological culture as described would certainly benefit from further field evaluation employing dedicated recto/vaginal swabs in order to fully establish diagnostic and analytical characteristics. Efforts are currently underway to address this.

Whilst this study utilised a GENIE® II fluorometer (cost £5000) to conduct isothermal amplification and specific detection of GBS in real time it would be entirely feasible to perform GBS LAMP testing using a simple heating block in conjunction with simple visual endpoint interpretation. This represents a commonly applied and valid approach for LAMP testing particularly in resource limited settings [27,38,39].

In conclusion the data presented demonstrates the efficacy of employing LAMP on minimally processed vaginal swabs for rapid direct molecular detection of GBS. The simplicity of the direct LAMP protocol combined with described sensitivity and specificity are desirable characteristics amenable to POC use in resource limited settings.

Acknowledgments

This project was funded by Northern Ireland Public Health Agency. RB was funded by an Academic Research Fellowship (Queen’s University Belfast). JPMcK, PVC, and DJF hold share options in HiberGene Diagnostics Ltd.

References

- 1. Blancas D, Santin M, Olmo M, Alcaide F, Carratala J et al.** Group B streptococcal disease in nonpregnant adults: incidence, clinical characteristics, and outcome. *Eur J Clin Microbiol Infect Dis* 2004;23:168–73.
- 2. Heath PT, Balfour G, Weisner AM, Efstratiou A, Lamagni TL et al.** Group B streptococcal disease in UK and Irish infants younger than 90 days. *Lancet* 2004;363:292-294.
- 3. Levent F, Baker CJ, Rench MA, Edwards MS.** Early Outcomes of Group B Streptococcal Meningitis in the 21st Century. *Pediatr Infect Dis J* 2010;29:1009–1012.
- 4. Verani, J.R., McGee, L. & Schrag, S.J. (2010).** Prevention of Perinatal Group B Streptococcal Disease: Revised Guidelines from CDC, 2010. *MMWR Recomm Rep* 2010;**59** (RR-10), 1-36.
- 5. Edmond KM, Kortsalioudaki C, Scott S, Schrag SJ, Zaidi AK et al.** Group B streptococcal disease in infants aged younger than 3 months: systematic review and meta-analysis. *Lancet* 2012;379:547-556.
- 6. Regan JA, Klebanoff MA, Nugent, R.P., Eschenbach, D.A., Blackwelder, W.C., Lou, Y., Gibbs, R.S., Rettig, P.J., Martin, D.H. & Edelman, R. (1996).** Colonization with group B streptococci in pregnancy and adverse outcome. *Am J Obstet Gynecol* 1996;**174**, 1354-1360.

- 374 **7. Goodman JR, Berg RL, Gribble RK, Meier PR. Fee RC *et al.*** Longitudinal Study of
375 Group B Streptococcus Carriage in Pregnancy. *Infect Dis Obstet Gynecol* 1997;5:237-243.
376
- 377 **8. Hansen SA, Uldbjerg N, Kilian M, Sørensen UBS.** Dynamics of *Streptococcus*
378 *agalactiae* Colonization in Women during and after Pregnancy and in Their Infants. *J Clin*
379 *Microbiol* 2004;42:83-89.
380
- 381 **9. Colbourn T, Gilbert R.** An overview of the natural history of early onset group B
382 Streptococcal disease in the UK. *Early Hum Dev* 2007;83:149–56.
383
- 384 **10. Gilbert R.** Prenatal screening for group B streptococcal infection: gaps in the evidence.
385 *Int J Epidemiol* 2004;33:2–8.
386
- 387 **11. Lin FYC, Weisman LE, Azimi P, Young AE, Chang K *et al.*** Assessment of
388 Intrapartum Antibiotic Prophylaxis for the Prevention of Early-onset Group B Streptococcal
389 Disease. *Pediatr Infect Dis J* 2011;**30**, 759-763.
390
- 391 **12. Yancey, M.K., Schuchat, A., Brown, L.K., Ventura, V.L. & Markenson, G.R. (1996).**
392 The accuracy of late antenatal screening cultures in predicting genital group B streptococcal
393 colonization at delivery. *Obstet Gynecol* 1996;**88**, 811-815.
394
- 395 **13. Davies HD, Miller MA, Faro S, Gregson D, Kehl SC *et al.*** Multicenter study of a rapid
396 molecular-based assay for the diagnosis of group B streptococcus colonization in pregnant
397 women. *Clin Infect Dis* 2004;39:1129-1135.

- 398 **14. Miller SA, Deak E, Humphries R.** Comparison of the AmpliVue®, BD MAX™
399 System, and illumigene® Molecular Assays for the Detection of Group B Streptococcus in
400 Antenatal Screening Specimens. *J Clin Microbiol* 2015;53:1938-1941.
- 401
- 402 **15. Poncelet-Jasserand E, Forges F, Varlet MN, Chauleur C, Seffert P et al.** Reduction of
403 the use of antimicrobial drugs following the rapid detection of *Streptococcus agalactiae* in
404 the vagina at delivery by real-time PCR assay. *BJOG* 2013;120:1098-1109.
- 405
- 406 **16. Gavino M, Wang E.** A comparison of a new rapid real-time polymerase chain reaction
407 system to traditional culture in determining group B streptococcus colonization. *Am J Obstet*
408 *Gynecol* 2007;97:388.e1-388.e4.
- 409
- 410 **17. Bergseng H, Bevanger L, Rygg M, Bergh K.** Real-time PCR targeting the sip gene for
411 detection of group B streptococcus colonization in pregnant women at delivery. *J Med*
412 *Microbiol* 2007;56:223-228.
- 413
- 414 **18. de Zoysa A, Edwards K, Gharbia S, Underwood A, Charlett A et al.** Non-culture
415 detection of *Streptococcus agalactiae* (Lancefield group B *Streptococcus*) in clinical samples
416 by real-time PCR. *J Med Microbiol* 2012;61:1086-1090.
- 417
- 418 **19. Park JS, Cho DH, Yang JH, Kim MY, Shin SM et al.** Usefulness of a Rapid Real-time
419 PCR Assay in Prenatal Screening for Group B Streptococcus Colonization. *Ann Lab Med*
420 2013;33:39-44.

- 421 **20. Couturier BA, Weight T, Elmer H, Schlager R.** Antepartum Screening for Group B
422 *Streptococcus* by Three FDA Cleared Molecular Tests and Effect of Shortened Enrichment
423 Culture on Molecular Detection Rates. *J Clin Microbiol* 2014;52:3429-3432.
- 424
- 425 **21. Francois P, Bento M, Hibbs J, Bonetti EJ, Boehme CC et al.** Robustness of Loop-
426 Mediated Isothermal Amplification Reaction for Diagnostic Applications. *FEMS Immunol*
427 *Med Microbiol* 2011;62:41-48.
- 428
- 429 **22. Niemz A, Ferguson TM, Boyle DS.** Point-of-care nucleic acid testing for infectious
430 diseases. *Trends Biotechnol* 2011;29:240-250.
- 431
- 432 **23. Heiniger EK, Buser JR, Mireles L, Zhang X, Ladd PD et al.** Comparison of point-of-
433 care-compatible lysis methods for bacteria and viruses. *J Microbiol Methods* 2016;128:80-87.
- 434
- 435 **24. Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K et al.** Loop
436 mediated isothermal amplification of DNA. *Nucleic Acids Res* 2000;28:e63.
- 437
- 438 **25. Nagamine K, Hase T, Notomi T.** Accelerated reaction by loop mediated isothermal
439 amplification using loop primers. *Mol Cell Probes* 2002;16: 223-229.
- 440
- 441 **26. Tomita, N., Mori, Y., Kanda, H. & Notomi, T. (2008).** Loop mediated isothermal
442 amplification (LAMP) of gene sequences and simple visual detection of products. *Nat*
443 *Protocols* 2008;3, 877-882.

27. Bourke TW, McKenna JP, Coyle PV, Shields MD, Fairley DJ. Diagnostic accuracy of loop-mediated isothermal amplification as a near-patient test for meningococcal disease in children: an observational cohort study. *Lancet Infect Dis* 2015;15:552-558.
28. Kaneko H, Kawana T, Fukushima E, Suzutani T. Tolerance of loop-mediated isothermal amplification to a culture medium and biological substances. *J Biochem Biophys Methods* 2007;70:499-501.
29. Kim DW, Kilgore PE, Kim EJ, Kim SA, Anh DD *et al.* The Enhanced Pneumococcal LAMP Assay: A Clinical Tool for the Diagnosis of Meningitis Due to *Streptococcus pneumoniae*. *PLoS ONE* 2012;7(8): e42954.
30. Huy NT, Hang TT, Boamah D, Lan NT, Van Thanh P *et al.* Development of a single-tube loop-mediated isothermal amplification assay for detection of four pathogens of bacterial meningitis. *FEMS Microbiol Lett* 2012;337:25-30.
31. Kimura K, Yanagisawa H, Wachino J, Shibayama K, Arakawa Y. Rapid and Reliable Loop-Mediated Isothermal Amplification Method for Detecting *Streptococcus agalactiae*. *Jpn J Infect Dis* 2013;66:546-548.
32. Suebsing, R., Kampeera, J., Tookdee, B., Withyachumnarnkul, B., Turner, W. & Kiatpathomchai, W. (2013). Evaluation of colorimetric loop-mediated isothermal amplification assay for visual detection of *Streptococcus agalactiae* and *Streptococcus iniae* in tilapia. *Lett Appl Microbiol* 2013;57, 317-324.

33. Slifkin, M. & Cumbie, R. (1987). Serogrouping Single Colonies of Beta-Hemolytic Streptococci with Achromopeptidase Extraction. *J Clin Microbiol* 1987;**25**, 1555-1556.

34. Leonard RB, Carroll KC. Rapid lysis of Gram positive cocci for Pulse field Gel Electrophoresis Using Achromopeptidase. *Diagn Mol Pathol* 1997;6:288-291.

35. Burke R, McKenna JP, Cox C, Coyle PV, Shields MD et al. A comparison of different pre-lysis methods and extraction kits for recovery of *Streptococcus agalacticae* (Lancefield group B Streptococcus) DNA from whole blood. *J Microbiol Methods* 2016;129:103-108.

36. Mahony J, Chong S, Bulir D, Ruyter A, Mwawasi K et al. Multiplex loop-mediated isothermal amplification (M-LAMP) assay for the detection of influenza A/H1, A/H3 and influenza B can provide a specimen-to-result diagnosis in 40 min with single genome copy sensitivity. *J Clin Virol* 2013;58:127-131.

37. Dunne WM Jr, Holland-Staley CA. Comparison of NNA agar culture and selective broth culture for detection of group B streptococcal colonisation in women. *J Clin Microbiol* 1998;36: 2298-2300.

38. McKenna JP, Fairley DJ, Shields MD, Cosby L, Wyatt DE et al. Development and Clinical Validation of a Loop mediated isothermal amplification (LAMP) method for the rapid detection of *Neisseria meningitidis*. *Diagn Microbiol Infect Dis* 2011;69:137-144.

491 **39. Hayashida K, Kajino K, Hachaambwa L, Namangala B, Sugimoto C.** Direct Blood
492 Dry LAMP: A Rapid, Stable, and Easy Diagnostic Tool for Human African Trypanosomiasis.
493 *PLoS Negl Trop Dis* 2015;9:e0003578. doi:10.1371/journal.pntd.0003578

494

Table 1 LAMP (Indirect & direct) and culture testing performance for clinical swabs (n=157) compared to reference GBS qPCR method. For LAMP the mean, maximum, minimum and median time to achieve LAMP positivity in minutes: seconds (T_T) and specific melt temperatures (T_m) in degrees centigrade are shown.

GBS qPCR	LAMP (INDIRECT)		LAMP (DIRECT)		Culture	
	Positive	Negative	Positive	Negative	Positive	Negative
Positive	43	0	41	2	36	7
Negative	0	114	0	114	0	114
Sensitivity	100%		95.4%		100%	
Specificity	100%		100%		94.2%	
PPV	100%		100%		83.7%	
NPV	100%		98.3%		100%	
(LR+)	Infinity		Infinity		Infinity	
(LR-)	0.00		0.05		0.16	
Mean T _T (T _m)	09.61 mins (86.66 ⁰ C)		14.69 mins (86.68 ⁰ C)		NA	
Max T _T (T _m)	29.00 mins (87.60 ⁰ C)		26.00 mins (87.60 ⁰ C)		NA	
Min T _T (T _m)	05.00 mins (85.86 ⁰ C)		08.15 mins (85.75 ⁰ C)		NA	
Median T _T (T _m)	08.30 mins (86.39 ⁰ C)		14.00 mins (86.40 ⁰ C)		NA	

NA; Not applicable

Table 2 Individual vaginal swab specimens identified either as GBS positive by culture, qPCR, indirect and/or direct LAMP methods. For qPCR genome copies per ml detected are shown. For LAMP the time to achieve positivity in minutes T_T and specific melt temperatures (T_m) in degrees centigrade are shown

Specimen Code	Culture Result	qPCR Result copies/ml	Indirect LAMP Result T _T mins (T _m)	Direct LAMP Result T _T mins (T _m)
334	+	+ 2.4E06	+ 08.20 (87.20 ⁰ C)	+ 17.50 (87.60 ⁰ C)
337	+	+ 5.6E03	+ 11.50 (87.60 ⁰ C)	+ 22.20 (87.60 ⁰ C)

354	+	+ 1.1E09	+ 05.30 (87.30 ⁰ C)	+ 09.51 (87.40 ⁰ C)
355	-	+ 4.2E03	+ 10.30 (87.11 ⁰ C)	+ 21.00 (87.10 ⁰ C)
370	+	+ 9.0E04	+ 09.30 (87.60 ⁰ C)	+ 19.00 (87.25 ⁰ C)
372	-	+ 6.3E03	+ 15.15 (87.41 ⁰ C)	+ 16.00 (86.00 ⁰ C)
373	+	+ 2.2E06	+ 08.00 (87.10 ⁰ C)	+ 14.50 (87.10 ⁰ C)
374	+	+ 2.5E09	+ 05.00 (87.30 ⁰ C)	+ 09.30 (87.41 ⁰ C)
376	+	+ 3.6E04	+ 09.30 (87.43 ⁰ C)	+ 20.15 (87.60 ⁰ C)
380	+	+ 2.8E04	+ 12.00 (87.09 ⁰ C)	+ 24.15 (87.20 ⁰ C)
382	+	+ 4.6E03	+ 17.15 (87.55 ⁰ C)	+ 21.00 (87.40 ⁰ C)
392	-	+ 5.5E03	+ 12.00 (87.40 ⁰ C)	+ 14.50 (86.61 ⁰ C)
395	+	+ 2.4E04	+ 11.45 (87.35 ⁰ C)	+ 20.30 (87.55 ⁰ C)
403	+	+ 1.1E04	+ 09.45 (86.04 ⁰ C)	+ 22.45 (87.38 ⁰ C)
405	-	+ 1.9E03	+ 12.15 (86.12 ⁰ C)	-
408	+	+ 1.7E05	+ 09.00 (86.09 ⁰ C)	+ 26.00 (87.44 ⁰ C)
425	+	+ 5.6E07	+ 07.00 (85.86 ⁰ C)	+ 08.45 (85.75 ⁰ C)
426	-	+ 3.7E04	+ 12.15 (86.14 ⁰ C)	+ 15.45 (86.12 ⁰ C)
427	+	+ 2.8E06	+ 08.30 (86.11 ⁰ C)	+ 11.30 (85.92 ⁰ C)
446	+	+ 1.5E08	+ 08.00 (87.00 ⁰ C)	+ 17.45 (87.11 ⁰ C)
448	-	+ 7.6E06	+ 09.45 (86.91 ⁰ C)	+ 20.00 (87.02 ⁰ C)
451	+	+ 1.3E08	+ 07.15 (87.00 ⁰ C)	+ 15.45 (87.21 ⁰ C)
453	-	+ 1.9E03	+ 29.00 (87.14 ⁰ C)	-
685	+	+ 2.2E08	+ 07.00 (86.30 ⁰ C)	+ 12.30 (86.50 ⁰ C)
687	+	+ 1.0E08	+ 07.20 (86.30 ⁰ C)	+ 11.00 (86.51 ⁰ C)
695	+	+ 9.0E04	+ 11.15 (86.49 ⁰ C)	+ 17.45 (86.39 ⁰ C)
711	+	+ 1.1E08	+ 06.50 (86.20 ⁰ C)	+ 10.30 (86.38 ⁰ C)
729	+	+ 5.6E07	+ 07.21 (86.20 ⁰ C)	+ 13.15 (86.43 ⁰ C)
732	+	+ 2.6E04	+ 11.50 (86.50 ⁰ C)	+ 18.15 (86.40 ⁰ C)
736	+	+ 9.4E05	+ 08.00 (86.30 ⁰ C)	+ 11.30 (86.30 ⁰ C)
755	+	+ 5.2E07	+ 07.00 (86.28 ⁰ C)	+ 10.45 (86.40 ⁰ C)
782	+	+ 2.4E06	+ 08.45 (86.49 ⁰ C)	+ 11.10 (86.30 ⁰ C)
789	+	+ 1.6E04	+ 17.00 (86.39 ⁰ C)	+ 19.00 (86.20 ⁰ C)
793	+	+ 2.3E07	+ 08.00 (86.49 ⁰ C)	+ 11.00 (86.40 ⁰ C)
825	+	+ 4.5E06	+ 08.15 (86.29 ⁰ C)	+ 10.45 (86.26 ⁰ C)
826	+	+ 1.1E07	+ 08.30 (86.29 ⁰ C)	+ 09.45 (86.30 ⁰ C)
831	+	+ 9.3E07	+ 06.45 (86.39 ⁰ C)	+ 09.00 (86.40 ⁰ C)

843	+	+ 4.2E08	+ 06.15 (86.29 ⁰ C)	+ 08.15 (86.36 ⁰ C)
864	+	+ 2.0E08	+ 06.30 (86.29 ⁰ C)	+ 08.15 (86.33 ⁰ C)
870	+	+ 7.0E05	+ 09.30 (86.39 ⁰ C)	+ 14.00 (85.80 ⁰ C)
872	+	+ 8.5E06	+ 08.00 (86.00 ⁰ C)	+ 10.15 (85.80 ⁰ C)
876	+	+ 5.2E07	+ 07.15 (86.29 ⁰ C)	+ 11.00 (86.42 ⁰ C)
891	+	+ 1.8E06	+ 08.30 (86.29 ⁰ C)	+ 11.00 (86.39 ⁰ C)
Totals	36	43	43	41

505